

Introducing the NEW
MONOJECT

Magellan*

Safety Needle



Discover a Whole
NEW World
of Needle Safety.



- **Easy one-handed safety activation**
 - Minimizes training and in-service
- **Needle-based engineering control**
 - Promotes clinical flexibility
- **Robust safety shield covers entire needle**
 - Maximizes protection from sharps injury
- **Integral safety shield**
 - Eliminates ability to remove safety device
- **Pre-attached syringe combinations**
 - Saves time and improves efficiency
- **Complete product line**
 - Simplifies the standardization to safety

tyco
Healthcare

Kendall

MONOJECT MAGELLAN* Safety Needle

Needle Only			Pre-Attached Combinations		
Code	Description	Packaging	Code	Description	Packaging
8881850810	18 x 1"	50/Bx, 500/Cs	8881833810	3cc 18 x 1"	50/Bx, 400/Cs
8881850815	18 x 1 1/2"	50/Bx, 500/Cs	8881833010	3cc 20 x 1"	50/Bx, 400/Cs
8881850910	19 x 1"	50/Bx, 500/Cs	8881833015	3cc 20 x 1 1/2"	50/Bx, 400/Cs
8881850915	19 x 1 1/2"	50/Bx, 500/Cs	8881833110	3cc 21 x 1"	50/Bx, 400/Cs
8881850010	20 x 1"	50/Bx, 500/Cs	8881833115	3cc 21 x 1 1/2"	50/Bx, 400/Cs
8881850015	20 x 1 1/2"	50/Bx, 500/Cs	8881833210	3cc 22 x 1"	50/Bx, 400/Cs
8881850158	21 x 5/8"	50/Bx, 500/Cs	8881833215	3cc 22 x 1 1/2"	50/Bx, 400/Cs
8881850110	21 x 1"	50/Bx, 500/Cs	8881833310	3cc 23 x 1"	50/Bx, 400/Cs
8881850115	21 x 1 1/2"	50/Bx, 500/Cs	8881833558	3cc 25 x 5/8"	50/Bx, 400/Cs
8881850210	22 x 1"	50/Bx, 500/Cs	8881833510	3cc 25 x 1"	50/Bx, 400/Cs
8881850215	22 x 1 1/2"	50/Bx, 500/Cs	8881866810	6cc 18 x 1 1/2"	50/Bx, 200/Cs
8881850358	23 x 5/8"	50/Bx, 500/Cs	8881866015	6cc 20 x 1 1/2"	50/Bx, 200/Cs
8881850310	23 x 1"	50/Bx, 500/Cs	8881866110	6cc 21 x 1"	50/Bx, 200/Cs
8881850558	25 x 5/8"	50/Bx, 500/Cs	8881866115	6cc 21 x 1 1/2"	50/Bx, 200/Cs
8881850510	25 x 1"	50/Bx, 500/Cs	8881866215	6cc 22 x 1 1/2"	50/Bx, 200/Cs
			8881822810	12cc 18 x 1"	50/Bx, 200/Cs
			8881822015	12cc 20 x 1 1/2"	50/Bx, 200/Cs
			8881822115	12cc 21 x 1 1/2"	50/Bx, 200/Cs
			8881822110	12cc 21 x 1"	50/Bx, 200/Cs
			8881822215	12cc 22 x 1 1/2"	50/Bx, 200/Cs



Sturdy sheath protects needle prior to use.



Activate MONOJECT MAGELLAN safety device using either thumb or finger.



Once locked, the safety shield completely encases the needle.



Activate MONOJECT MAGELLAN safety device on a hard surface.

*MONOJECT MAGELLAN is a trademark of Sherwood Services AG.

tyco
Healthcare

Kendall

15 Hampshire Street
Mansfield, MA 02048
1-800-962-9888
508-261-8000

www.kendallhq.com

For International locations please see

www.tycohealthcare.com

H-5350 5M 4/03